**TPS DS**

**TPS**

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### THINH PHAT SCALES JONT STOCK COMPANY

Address : 57 Street D1, Ward 25, Binh Thanh District  
HoChi Minh City Viet Nam

Tel : (08) 62.999.111 , (08) 62.888.666 , (08) 3512.7131(12 line)

Fax : (08) 3512.7699 , Order (08) 62.999.666

Website : [www.canthinhphat.com](http://www.canthinhphat.com) - [www.canvina.com](http://www.canvina.com)

Email : [hcm@canvina.com](mailto:hcm@canvina.com) - [vo@canvina.com](mailto:vo@canvina.com)

**MAIN FEATURES AND SPECIFICATION****CAUTION**

- 1.1 Please turn off the charge mode if do not use the hydrogens battery of nickel or Lead — aid battery.
- 1.2 Avoid installing the balance in direct sunlight-this may cause poor performance or malfunctions.
- 1.3 Do not mix different type of batteries. Replace all the Batteries at the same time.
- 1.4 If the balance is not to be used for a long period of time, remove all batteries from the battery compartment to avoid leakage which may cause damage to the instrument.
- 1.5 Avoid overloading or dropping material onto the platform as that could damage the balance.
- 1.6 Do not spill liquids on the balance as it is not water-resistant. Liquids may damage the case and if it gets inside the balance it may cause damage to the electronics.
- 1.7 Material that has a static electric charge could influence the weighing. Discharge the static electricity of the samples, if possible. Another solution to the problem is to wipe both sides of the pan and the top of the case with an antistatic agent.

**PREPARATION**

1. Place your balance on a firm and level table and then level the balance by adjusting the feet to center the bubble in the spirit Level.
2. Avoid placing the balance in direct strong sunlight or draught that may cause poor performance or malfunctions.
3. Use the individual socket to avoid the interfere of other wirings.
4. Turn on the balance while there is no load on the pan.
5. Place the weight on the center of the pan for accurate weighing. Make sure the weight do not over the Max capacity as it could damage the load cell inside.
6. We suggest to warm — up the balance for 5 minutes for accurate weighing.
7. Calibration may be required before weighing or when the balance is initially installed or moved from a location.
8. Please replace the battery when the low battery  mark comes.

Introduction : This is a multi-function weighing balance that offer flexibility together with good features!

Main features:

- A: Automatic zero point tracking function.
- B: With multi-rang calibration.
- C: Tare function.
- D: Many units can be selection.

E: Simple counting function and percentage function are applicable. F: Fast, stable and accurate display reading with simple operation.

**PRECAUTIONS:**

1. Full charge the battery after unpacking the scale.  
 Recharge the battery: *When battery symbol appears on the LCD display, charge the battery with AC adapter, the balance of charge will light up in red, when it becomes green means charge completed. (It takes about 8 hours to full charge the battery.)*
2. Install the equipment on a level and stable surface.
3. Do not install the equipment near the air conditioning or a vibrating machine.
4. Install the equipment in an environment with steady temperature(0°C ~40°C), prevent from rapid temperature changes.
5. Independent AC adapter for this equipment is recommended, check the
6. Warm up the equipment for 15 minutes before use.
7. Avoid overloading or dropping material onto the platform.

**SPECIFICATIONS for size 120mm**


Model	DS303	DS603	DS302	DS602	DS1002	DS1502
Maximum Capacity	300g	600g	300g	600g	1000g	1500g
Division (d)	0.001g	0.001g	0.01g	0.01g	0.01g	0.01g
Readability (e)	0.05 g	0.05 g	0.1 g	0.1 g	0.1 g	0.1 g
Min	0.01g	0.01g	0.2g	0.2g	0.2g	0.2g
Display Resolution	1/60000					
Response	3~5 seconds					
Display	LCD display with 6 digits / white backlight Model: TPS Seri C_DOUBLE_HT1621LCD DB05053A Revision: 1.1 produce by TPS Corporation					
Loadcell	Model: LAB-B (B is capacity of loadcell) produce by TPS Corporation					
Unit	g, ct, lb, oz, tlt (can be set by soft)					
Overload	150 %					
Appearance	ABS, with stainless steel pan, breeze shield					
Calibration weight	3g – 1500g					
Power source	Battery –DC6V 1.3Ah Adaptor –input: AC 110-240V output: 12V/500mA					
Temperature	preservation: -10°C~+50 °C; Work: 0°C~+40 °C					
Moisture	preservation 5%R.H~90%R.H; Work: 10%R.H~80%R.H					
Dimensions (LxDxH)	(182 x 283 x 70) mm					
Pan-size (LxD)	120mm					
Net Weight	1.5 kg					

Connect the balance with PC or printer, this function is effective.  
 Under weighing mode, press the **[PRINT]** key to print out percent weight value.

**SPECIFICATIONS for size (160x160)mm**


Model	DS2002	DS3002	DS6002	DS3001	DS6001	DS7501
Maximum Capacity	2kg	3kg	6kg	3kg	6kg	7.5kg
Division (d)	0.01g	0.01g	0.01g	0.1g	0.1g	0.1g
Readability (e)	0.5 g	0.5 g	0.5 g	0.5 g	0.5 g	0.5g
Min	0.2g	0.2g	0.2g	2g	2g	2g
Display Resolution	1/60000					
Response	3~5 seconds					
Display	LCD display with 6 digits / white backlight Model: TPS Seri C_DOUBLE_HT1621LCD DB05053A Revision: 1.1 produce by TPS Corporation					
Loadcell	Model: LAB-B (B is capacity of loadcell) produce by TPS Corporation					
Unit	g, ct, lb, oz, tlt (can be set by soft)					
Overload	150 %					
Appearance	ABS, with stainless steel pan, breeze shield					
Calibration weight	3g – 1500g					
Power source	Battery –DC6V 1.3Ah Adaptor –input: AC 110-240V output: 12V/500mA					
Temperature	preservation: -10°C~+50 °C; Work: 0°C~+40 °C					
Moisture	preservation 5%R.H~90%R.H; Work: 10%R.H~80%R.H					
Dimensions (LxDxH)	(182 x 283 x 70) mm					
Pan-size (LxD)	(160 x 160)mm					
Net Weight	1.7 kg					

Connect the balance with PC or printer, this function is effective.  
 Under weighing mode, press the **[PRINT]** key to print out percent weight value.

**DISPLAY AND KEY DESCRIPTION**

**LCD display**

Indication sign: stable zero Tare low battery

**Start to use: Accessories:**

1. adapter
2. balance body
3. plastic pan
4. SS pan
5. wind shield

**Turn on the indicator**

- Press the ON/OFF] key.
- Wait till the screen displays “0.00”. If it doesn’t shows “0.00”, press [ZERO] key to enter weighing mode.

**KEYPAD AND PANEL INTRODUCTION**

	Lock the display weight value		Change Unit g, ct, lb, oz, tlt
	Change Weight/PCS/% functions		TARE
	Print out the data		Reset to zero

**FUNCTIONS AND OPERATION**

**1. Simple weighing**

1. Press the on/off] key or Turn on the power switch on the right side of balance.
2. The balance power on with beep and backlight flash.
3. Wait till the screen displays “0.00g”.Press [ZERO] key to enter weighing mode if the screen doesn’t display “0.00g”.
4. Put the object on weighing pan then read the weight within 2~3 seconds.

**2.Tare function:**

- Put the container (object) on weighing pan, press [TARE] to deduct the container(object) weight, the screen displays “0.00g”.
- Put more objects on the weighing pan, the screen displays weight of new objects excluding weight of container.

**3.Units shift and setting**

Under weighing mode, press [UNITS] key to change to the unit you need. Open more units or close any units, please contact your local agent.

#### 4. Simple counting setting

1. Under weighing mode, press [PCS] choose the simple counting function. Put on the weights, if no unit weight setting before, screen will display “0 pcs”.
2. Press [UNIT] key, enter quantity setting, The screen displays “N-XXX” (“XXX” stands for sample quantity.)
3. Press [UNIT] key, move the digit.
4. Press [TARE] key, add up the percent digits.
5. Press [ZERO] key, reduce the percent digits.
6. Press [PCS] key, save the above setting, you can begin to use this function.

#### 5. Simple percent counting

In counting mode, press [MODE] to change % mode or weighting mode.

- Under weighing mode, press [MODE], choose the % function. Put on the weights, screen will display “S-100 %”.
- Press [MODE] key, save the above setting, you can begin to use this function.

\*\* LED model do not have % mode, so only change to weighting mode.

#### 6. Hold function

1. Press once to lock the percent weight value.
2. Press again to cancel lock status.

\*\* This function usually set to close, if you want to use this function, please contact your local agent.

#### 7. Printing function

connect the balance with PC or printer, this function is effective.

Under weighing mode, press the [PRINT] key to print out percent weight value.

#### 8. Simple calibration

01	 <p>Press hold key 2 seconds to enter simple calibration</p>	<p><b>CAL</b></p> <p>Press [UNIT] to move the digits                  Press [TARE] key, add up the percent digits                  Press [ZERO] key, reduce the percent digits                  Put on the right weights (the same as the display weight).</p>
02	 <p>Press key</p>	<p>Wait the stable symbol appears for 2 seconds the buzzer will send out a short beep, that means simple calibration successfully.</p>

**ERROR MESSAGE AND TROUBLESHOOTING**

Error message

Message	Problem
ADO---	Over A/D resolution range
ADL---	Low A/D resolution range
OVR---	Over load (max: capacity +9e)
Battery symbol	Low battery alarm
Lo-bAt	Low battery, Need charge

Trouble shooting

Where	Possible Error message	Trouble shooting.
Power on	ADO---	Check and remove the object from weighing pan or send back for
Power on	ADL---	send back for repair.
Power on	Battery symbol	Charge the battery.
Normal weighing mode	OVR---	Check if weighing object over the capacity +9e.
Normal weighing mode	Lo-bAt	Charge the battery.

---

**Hotline Call to our company TPS Corporation in Viet Nam**

Call	Tel
Hotline	(+84.8)62.888.666 / (+84.8) 62.999.111
Mr Cuong (24/24)	0084 915.999.111
Mr Xa	0084 974.000.333
Ms Cuong	0084 908.444.000
Delivery	0084 919.555.444

**THINH PHAT SCALES JONT STOCK COMPANY**

Address : 57 Street D1, Ward 25, Binh Thanh District  
HoChi Minh City Viet Nam

Tel : (08) 62.999.111 , (08) 62.888.666 , (08) 3512.7131 (12 line)

Fax : (08) 3512.7699 , Order (08) 62.999.666

Website : [www.canthinhphat.com](http://www.canthinhphat.com) - [www.canvina.com](http://www.canvina.com)

Email : [hcm@canvina.com](mailto:hcm@canvina.com) - [vo@canvina.com](mailto:vo@canvina.com)